

## Model : AD7013MT Dino-Lite Premier



- 5MP sensor
- Detachable interchangeable front caps
- Aluminum alloy case
- MicroTouch II
- Requires Motion Jpeg (MJPEG) codec to view at 5 megapixels (2592x1944)
- OS: Windows XP, Vista, 7
- Unit Weight: 155g
- Units Dimension 10cm (H) x 3.2cm (D)

### Product Description

Product Number	AD7013MT	Magnification	up to 200x
Warranty Period	1	Pixels	5M
Condition	New	Resolution	2592x1944
Packaging	Retail Box	Data Output	USB 2.0
Sensor	Enhanced Color CMOS	Video frame Rate	up to 30 frame-per-second
Illumination	Buit-in 8 white LEDs	Save File format	BMP, JPG, AVI
LED On/Off Function	Yes	Body length	4 inches / 10.2 cm
Body Diameter	1.26 inches / 3.2cm max.	Cable length	6 feet / 180 cm
Unit weight	155g		

### Product Applications

The Dino-Lite Premier AD7013MT digital handheld microscope has a 5 megapixel sensor allowing you to see and record every detail providing you an experience you never had with a digital handheld microscope. It also includes a valuable capability to remove the front cap and interchange it with another that fits the situation, such as the included open cap for general purposes and closed cap for extra protection against dust and debris.

The AD7013MT is appropriately incased with an aluminum alloy casing providing robust protection only fit to hold 5 megapixels of viewing power. It can magnify up to 200x magnification depending of working distance and includes our convenient MicroTouch II sensor for on microscope picture taking. It is also illuminated by built in bright white LED's to instantly illuminate objects.

The Dino-Lite digital handheld microscope experience is not complete without the included DinoCapture 2.0 software that allows you to take pictures, record video, and also annotate on images. This model also includes the measurement and calibration capabilities only available for professional models such as this. The AD7013MT with 5 megapixel sensor allows you to see fine detail in any magnification for an experience you never had with a digital handheld microscope of this kind.

Included Caps:  
FC-OC1 (open cap), FC-T-CC1 (closed cap)



5.0 megapixel sensor



Microtouch II



USB

**Metal**



White Light



Measurement Calibration