


Model : AM311H Dino-Lite Earscope

- Health Care Digital Microscope
- 20x~30x
- Built-in 4 white LEDs
- 640x480 pixels. (VGA)
- USB2.0
- Windows Vista, XP, 2000, or Windows Server 2003
- MAC OS (Microtouch function no supported)



Product Description

Product Number	AM311H	Magnification	20x~30x
Warranty Period	1	Pixels	0.3M
Condition	New	Resolution	640x480
Packaging	Retail Box	Data Output	USB 2.0
Sensor	Enhanced Color CMOS	Video frame Rate	up to 30 frame-per-second
Illumination	built-in 4 white LED's	Save File format	BMP, JPG, AVI
Body length	4 inches / 10.2 cm	Body Diameter	1.26 inches / 3.2cm max.
Cable length	6 feet / 180 cm	Total package weight	280g +/- 20g

Product Applications

The Dino-Lite digital handheld ear scope is a portable, easy to use, and yet light weight ear scope that fits on the palm of your hands. The AM311H provides users the ability to view and record images of the ear canal/ eardrums with ease.

Included with each Dino-Lite handheld digital microscope is a copy of the "DinoCapture" software, enabling you to "capture" either a photograph, video, or time-lapsed video within a few simple clicks