



Model : AM311ST Dino-Lite

- Handheld Digital Microscope
- 0.3M / Resolution 640x480
- USB 2.0 Output
- 10~200X Continuous Magnification
- LED on/off controlled by software
- Built in Micro Touch (Trigger Button)

Microtouch

Product Description

Product Number	AM311ST Dino-Lite	Magnification	10~200X Variable Magnification
Warranty Period	1	Pixels	0.3M
Condition	New	Resolution	640x480
Packaging	Retail Box	Data Output	USB 2.0
Sensor	1/4" Color CMOS	Video frame Rate	up to 30 frames per second
Illumination	built-in 8 white LEDs	Save File format	BMP, JPG, AVI
Body length	4 inches / 10.2 cm	Body Diameter	1.26 inches / 3.2cm max.
Cable length	6 feet / 180 cm	Unit weight	90(g)



Plant



Plant 50X




Plant 200X



Printing Material



Printing 50X



Printing 200X

Product Applications

An amazing, easy to use handheld Micro Touch (Trigger Button) digital microscope with USB output interface. High resolution magnification at your fingertips! Detailed examination of any area you choose. The AM311ST handheld digital microscope provides sharp image quality and includes 8 built-in white LEDs with other wonderful yet functional features. This handheld digital microscope supports Windows XP, Windows Vista, Macintosh OS Tiger, and Macintosh OS Leopard. It is perfect for the home user who demands a simple USB handheld microscope to examine objects and take quality pictures. The AM311ST surely would give you high magnification and detailed images right from your fingertips.

The AM311ST is ideal for a wide variety of applications; such as in the classroom, workplace, and even at home. The AM311ST can be used for science, engineering, work, study, dermatology (skin/scalp), exam, detailed repair, assembly, quality control (i.e. electronics, mechanical, etc), hobbies (coin, stamp, watch, jewelry, etc), law enforcement (counterfeit ID, crime lab, etc), entertainment, reading aid, and much more. Just by using your imagination, the possibilities will be endless.