



## Model : AM4016T Dino-Lite Premier VGA Dino-Lite for VGA Monitor

- Resolution: 640x480pixel
- Microtouch for LED on/off
- Output:D-SUB
- Power Adaptor: 5V

Product Description			
Product Number	AM4016T	Magnification	200x
Warranty Period	1	Pixels	800x600pixels
Condition	New	Resolution	800x600pixels
Packaging	Retail Box	Data Output	D-SUB
Sensor	Enhanced Color CMOS	Video frame Rate	up to 30 frame-per-second
Illumination	Buit-in 8 white LEDs	Save File format	BMP, JPG, AVI
LED On/Off Function	Yes	Body length	4 inches / 10.2 cm
Body Diameter	1.26 inches / 3.2cm max.	Cable length	6 feet / 180 cm
Total package weight	280g +/- 20g		

### Product Applications

Enjoy direct connection to your computer monitor with the AM4016T Dino-Lite Premier VGA microscope that can connect directly to your computer monitor via VGA (D-SUB) connection. It can view 60fps at 800x600 resolutions and has enhanced low light visibility so you can see clearly in any situation.

With our unique patented design, the AM4016T can magnify objects up to 200x depending of working distance without changing a single lens. It also has a MicroTouch II sensor on the microscope so you can turn the built-in lights on or off. The AM4016T provides a versatile handheld microscope for your microscopy needs without a computer or a distracting television monitor. All you need other than the AM4016T is a computer monitor and your ready to go, providing low total system cost and focused use of our handheld microscope.

