



Model : AM4113RUT Dino-Lite Premier

- Iriscope
- 10X ~ 20X
- Enhanced 1.3 Megapixels
- Scroll lock
- MicroTouch II

Product Description

Product Number	AM4113RUT	Magnification	10X ~ 20X
Warranty Period	1	Pixels	1.3 Megapixels
Condition	New	Resolution	1280×1024 pixels. (VGA)
Packaging	Retail Box	Data Output	USB 2.0
Sensor	Enhanced Color CMOS	Video frame Rate	up to 30 frame-per-second
Illumination	8 white LED	Save File format	P, JPG, AVI
LED On/Off Function	Yes	Calibration Function	Yes
Body length	4 inches / 10.2 cm	Body Diameter	1.26 inches / 3.2cm max.
Cable length	6 feet / 180 cm	Unit	90(g)



Iris



1.3 megapixel sensor

Crystal Clear Images
1.3 Megapixels

Product Applications



The AM4113RUT digital handheld Iriscope is a simple, lightweight, and portable tool to take pictures of the iris. Simply place the Iriscope to the eye and take a picture. The AM4113RUT views with a 1.3 megapixel sensor which lets you view and record images with great detail. The AM4113RUT digital Iriscope is illuminated by two off set LED's to reduce glare.

The AM4113RUT digital handheld Iriscope comes with the DinoCapture software. With this innovative software you can view, take pictures, and annotate images right on your computer.



For Windows users the AM4113RUT handheld digital microscope is featured with the unique MicroTouch snapshot feature. The MicroTouch (touch sensitive static sensor) feature on the handheld digital microscope gives users the option of taking a picture right from the scope itself. The MicroTouch was designed to minimize image blur while taking a picture with the microscope. Great for those moments when you need to take a picture, but having a hard time reaching the computer.

This unit provides high resolution, professional tools, ease of use, and flexibility for the user.