




Model : AM423 Dino-Eye

- Eyepiece Digital Microscope
- AV /TV Output
- 1.3M / Resolution : 1280 x 1024
- Sensor: 1/3" Color CMOS
- Body diameter: 23 mm
- Body length: 52mm
- Aperature diameter: 23 mm



• AM423



• AM423



• AM423 with eyepiece



• Education/Science Lab



• Printing



• Eye rabbit testis

Product Description

Product Number	AM423	Magnification	none
Warranty Period	1	Pixels	1.3M
Condition	New	Resolution	1280x 1024
Packaging	Retail Box	Data Output	USB 2.0
Save File format	BMP, JPG, AVI	Video frame	Rate up to 30 frame-per-second
Sensor	1/3" Color CMOS	Aperature diameter	23 mm
Body diameter	23 mm	Body length	52 mm

Product Applications

The AM423 DinoEye digital eyepiece can help users upgrade their existing microscope, by replacing the eyepiece of a traditional microscope eyepiece with the DinoEye digital eyepiece.

The AM423 digital eyepiece uses an USB 2.0 interface that fits eyepiece slots with an inner diameter of 23mm.

You can use the AM423 digital eyepiece as a teaching tool, for digital imaging documentation, or simply for comfort on those long hours viewing into a microscope.

Included with each DinoEye digital eyepiece is a copy of the "DinoCapture" software, enabling you to "capture" either a photograph, video, or time-lapsed video within a few simple clicks.