



## Model : AM451 Dinolite Tabletop Digital Microscope

- Portable tabletop Digital Microscope
- 1.3M / Resolution 1280x1024
- USB(2.0) Output
- Three objective lens' 50x, 100x~200x, 400x~600x



### Product Description

Product Number	AM451	Magnification	50x 100x~200x / 400x~600x
Warranty Period	1	Pixels	1.3M
Condition	New	Resolution	1280 x 1024 pixel
Packaging	Retail Box	Data Output	USB 2.0
Sensor	Enhanced Color CMOS	Illumination	Oblique and Backside LED lighting
Cable length	6 feet / 180 cm	Total package weight	1,135g

### Product Applications

An amazing, easy to use portable desktop, three objective lens digital microscope with USB 2.0 output interface. High resolution magnification at your fingertips! Detailed examination of any area you choose. Model's AM451 is self illuminating with 4 LED's. - Ideal for a wide variety of practical applications in the classroom, workplace and even the home. Can be used for science and Laboratory Research, dermatology (skin/scalp) exam, Medical analysis, Education, Classroom, Teacher Aid as Project, Presentation, and more...

Included with each Dino-Lite digital microscope is a copy of the "DinoCapture" software, enabling you to "capture" either a photograph, video or time-lapsed video with a simple click!