

Model : AM7013MZT4 Dino-Lite Premier



- Handheld Microscope
- 5 Megapixel sensor
- Polarizer
- Aluminum alloy case
- MicroTouch II
- Unit Weight: 155g
- Units Dimension 10cm (H) x 3.2cm (D)
- MJPEG codec is required to run at 2592 x 1944 resolutions

Handheld Digital Microscope(USB)
400x~470x
Adjustable polarization
5MP sensor
With measurement and calibration
Aluminum alloy case

Product Description

Product Number	AM7013MZT	Magnification	up to 200x
Warranty Period	1	Pixels	5M
Condition	New	Resolution	2592x1944
Packaging	Retail Box	Data Output	USB 2.0
Sensor	Enhanced Color CMOS	Video frame Rate	up to 30 frame-per-second
Illumination	Buit-in 8 white LEDs	Save File format	BMP, JPG, AVI
LED On/Off Function	Yes	Body length	4 inches / 10.2 cm
Body Diameter	1.26 inches / 3.2cm max.	Cable length	6 feet / 180 cm
Unit weight	155g		

Product Applications

The Dino-Lite Premier AM7013MZT digital handheld microscope has a 5 megapixel sensor with adjustable polarizer abilities. Allowing you to observe things in high resolution and see things you couldn't with direct lighting. With up to 200x magnifications depending on working distance, polarization, and 2592x1944 of resolution; you will not miss detail under any magnification.

The AM7013MZT is appropriately incased with an aluminum alloy providing robust protection only fit to hold 5 megapixels of viewing power. It also includes our convenient Microtouch II sensor for on microscope picture taking and built in bright white LED's to instantly illuminate objects. The Dino-Lite digital handheld microscope experience is not complete without the included DinoCapture 2.0 software that allows you to take pictures, record video, and also annotate on images.

This model also includes the measurement and calibration capabilities only available for professional models such this. The AM7013MZT with 5 megapixels sensor and adjustable polarizer allows you to observe in overwhelming detail under any magnification and see more with its adjustable polarizer.

Notes:

- 5MP sensor
- Polarizer
- Aluminum alloy case
- MicroTouch II
- MJPEG codec is required to run at 2592 x 1944 resolutions
- Unit Weight: 155g
- Units Dimension : 10cm (H) x 3.2cm (D)
- Main USB is enough for most users but if experiencing problems plug in both USB



Metal

