



Model : AM7023 Dino-Eye

- Eyepiece Digital Microscope
- USB 2.0 Output
- 5M / Resolution 2592 x 1944
- Body diameter: 23 mm
- Body length: 84 mm
- Aperature diameter: 23 mm



Product Description

Product Number	AM7023	Magnification	none
Warranty Period	1	Pixels	5M pixels
Condition	New	Resolution	2592 x 1944
Packaging	Retail Box	Data Output	USB 2.0
Save File format	BMP, JPG, AVI	Video frame	Rate up to 30 frame-per-second
Sensor	1/3" Color CMOS	Aperature diameter	23 mm
Body diameter	23 mm	Body length	84 mm


Product Applications

The AM7023 DinoEye digital eyepiece can help users upgrade their existing microscope, by replacing the eyepiece of a traditional microscope eyepiece with the DinoEye digital eyepiece.

The AM7023 digital eyepiece uses an USB 2.0 interface that fits eyepiece slots with an inner diameter of 23mm.

You can use the AM7023 digital eyepiece as a teaching tool, for digital imaging documentation, or simply for comfort on those long hours viewing into a microscope.

Included with each DinoEye digital eyepiece is a copy of the "DinoCapture" software, enabling you to "capture" either a photograph, video, or time-lapsed video within a few simple clicks.




Malaria Parasites



Micro-Parasites



Mosquito



Bee