




Model : AMH-N5UT Dino-Lite Nailscope

- Health Care Digital Microscope
- 500x
- Built-in 8 white LEDs
- 1.3M pixels. (SXGA)
- USB2.0
- Windows Vista, XP, 2000, or Windows Server 2003
MAC OS (Microtouch function no supported)



Capillary



Capillary with
resolution oils



Capillary

Product Description

Product Number	AMH-N5UT	Magnification	500X
Warranty Period	1	Pixels	1.3M
Condition	New	Resolution	1280x1024 pixel
Packaging	Retail Box	Data Output	USB 2.0
Sensor	Enhanced Color CMOS	Illumination	built-in 8 white LED's
Body length	4 inches / 10.2 cm	Body Diameter	1.26 inches / 3.2cm max.
Cable length	6 feet / 180 cm	Total package weight	280g +/- 20g

Product Applications

The AMH-N5UT digital handheld nail microcirculation scope can be used for diagnosis of various kinds of diseases by checking the capillary vessel of your fingers. It is an ideal tool for health inspection and consultation. The digital handheld microscope views the capillaries at 500X. Images recorded and viewed are in 1.3 megapixels, which delivers clear and crisp images of the capillaries right onto your computer monitor. The AMH-N5UT comes with the DinoCapture software which lets you take pictures, view images, record video, and annotate images taken right from the AMH-N5UT onto your computer monitor.