





Model : AMH-RUT Dino-Lite Digital Iriscope


- Handheld Iris Digital Microscope
- 1.3M / Resolution 1280x1024
- 10~20X Continuous Magnification
- Achieves effective magnification values of up to 30x w/ digital zoom (only available in VGA mode - 640*480)
- USB 2.0 Output
- 2 LED on/off controlled by software
- w/Measurement



Iris Example 1



Iris Example 2



Iris Example 3

Product Description

Product Number	AMH-RUT	Magnification	10X ~ 20X
Warranty Period	1	Pixels	1.3M
Condition	New	Resolution	1280 x 1024 pixel
Packaging	Retail Box	Data Output	USB 2.0
Sensor	Enhanced Color CMOS	Illumination	built-in 8 white LED's
Body length	4 inches / 10.2 cm	Body Diameter	1.26 inches / 3.2cm max.
Cable length	6 feet / 180 cm	Total package weight	280g +/- 20g

Product Applications

The AMH-RUT digital handheld Iriscope is a simple, lightweight, and portable tool to take pictures of the iris. Simply place the Iriscope to the eye and take a picture. The AMH-RUT views with a 1.3 megapixel sensor which lets you view and record images with great detail. The AMH-RUT digital Iriscope is illuminated by two off set LED's to reduce glare.

The AMH-RUT digital handheld Iriscope comes with the DinoCapture software. With this innovative software you can view, take pictures, and annotate images right on your computer.