



Visioscope® BW 30

CK-Visioscope® BW 30

Images taken with a special UV-light source



These unique video cameras deliver very impressive images of skin, hair and scalp.

The cameras show a skin area of 6 x 8 mm and monitor:



Technical Data:

Dimensions: approx. 118 x 56 x 48 cm, weight: 250 g, Image size: 6x 8 mm; Power Supply: BNC plug: optional external power supply: input 100-240 V AC, 50/60 Hz, 0.3 A, output: 7.5 V DC, max. 1.5 A

Technical changes may be made without prior notice.