

Visioscope® Color 32

- skin texture (smoothness, wrinkles)
- desquamation (scaliness)
- skin impurities (reddening, pigmentation spots, acne, comedones etc.)
- hair structure, dandruff and condition of the scalp



With the special foils Sebufix®/Corneofix® also sebum and hydration levels of the skin can be shown.



CK-Visioscope® Color 32

Coloured images taken with white light.



The cameras must be connected to a video monitor. For the Visioscope® BW 30 there is also a Visioscope® Monitor Center available.

Technical Data:

Dimensions: approx. 118 x 56 x 48 mm, weight: 250 g, Image size: 6x 8 mm; Power Supply: BNC plug: optional external power supply: input 100-240 V AC, 50/60 Hz, 0.3 A, output: 7.5 V DC, max. 1.5 A; USB plug: via the computer

Technical changes may be made without prior notice.