

Multi Skin Test Center® MC 750

Due to new CK probe technology the device allows operation of all CK probes at the same time according to the demands of each user. The very comfortable and clear structured software simplifies the work with the system and allows to give individual product recommendations for skin care, hair care and sun screen products.

CK-measuring probes for the following parameters can be connected: moisture, oil, pH, pigmentation & redness, transepidermal waterloss (TEWL – the barrier function of the skin), skin temperature and elasticity.

In addition a special software for the use of dermatologists in order to save patient data and measuring values can be used with the same device. Individual treatments, evaluation and recommendations can be printed out. The software Dermacheck also allows to connect the Visioscope Color 32 camera via USB port to your computer to look at skin images and to save them.

Technical Data:

*Power supply: 100-240 V AC, 0.3 A, 50-60 Hz; Dimensions: 33 x 10 x 25 cm; Weight without probes: 3.2 kg
Device must be operated by an IBM-compatible computer
Technical changes may be made without prior notice.*

multitude of parameters possible

