



### Model : MS-09B RACK Mobile Stability Stand

- Stability meets portability
- 2 different bases: traditional flat and horseshoe base
- Provides stable viewing on the go
- Holster compatible with Dino-Lite handheld microscopes



#### Product Description

Product Number	MS-09B	Dimension	13.5cm (H)
Warranty Period	1	Type	Ultra-Portable Mini Stand
Condition	New	Unit Weight	186g

#### Product Applications

The MS-09B brings you a similar table top type of design as the MS35B stand with the portability of mobile stand. Because of its miniature size this stand will take up much less space when you need to minimize space in your office, workplace or demonstration environment.

The MS-09B stand comes complete with 2 different bases for different needs. First is a regular flat base with a slight 3x3 grid grooves to help center objects. Second is a 'horseshoe' type of base providing a direct view through the stand for various applications. Both bases have rubber grip pegs to help hold the stand in place and simultaneously not damage the surface you are placing the stand on. The MS-09B is another versatile tool in the expanding line of Dino-Lite accessories to help you maximize your potential.

Package includes:

- Holster
- Horseshoe base
- Flat table top base
- Mounting pole

