



Model : MS34B-R2 Table Top Stand

- Compact fine adjustment stand
- Lightweight and portable
- Provides stable viewing
- Integrated holster compatible with Dino-Lite handheld microscopes



Product Description

Product Number	MS34B-R2	Dimension	11cm (L) x 8cm (W) x 13.5cm (H)
Warranty Period	1	Type	Table top Stand
Condition	New	Unit Weight	262 g

Product Applications

The MS34B-R2 allows precision focus and stable viewing in a compact design.

With its precision vertical movement of the holster, the MS34B allows steady viewing while under higher magnifications. With its steady viewing, it also allows ease of use when fine focusing with the stand. Works exceptionally well with our Dino-Lite 500 x models like the AM413T5 for quick or prolonged and stationary use.

With its compact design, it will take up minimal space at home or at work and can be easily carried around with you. The MS34B-R2 is a great stand to assist focusing at high magnification and stable viewing wherever you a

