


Model : MS52B Table Top Clamp-on version

- Flexible articulating desktop stand
- Horizontal & Vertical dist: 23~36 inches (621~921mm)
- Top swivel: 360 degrees
- Weight: 4 lb



Product Description

Product Number	MS52B	Dimension	14" x 6 1/2" x 2"
Warranty Period	1	Type	Table Top Clamp On
Condition	New		

Product Applications

The MS52B is a jointed articulating flex-arm type stand with holster for use with the Dino-Lite hand held digital microscopes. The precision machine crafted components provide smooth adjustment and unlimited horizontal positioning in a well balanced ergonomic design. The MS52B can be easily clamped onto and secured to any workstation or table top. Manufactured from durable, lightweight recycled aluminum and able to be folded flat makes the versatile MS52B stand conveniently portable.

Note: The available MS52BA2 option includes the MS52B flex-arm stand with an additional rigid vertical floor stand and base assembly for use when a workstation table top is not convenient or available.

Table Top Clamp-on Rack Flexible articulating desktop rack Horizontal & Vertical Dist: 621~921mm Top swivel: 360 degree