



Model : MSAA502 Multi-purpose Adapter

- Multi-purpose Adapter
- 124mm (h) x 84 mm (d)



Product Description			
Product Number	MSAA502	Dimension	12.5cm (H) x 8.5cm (D)
Warranty Period	1	Type	Multi-purpose Adapter
Condition	New	Unit Weight	72g

Product Applications

The MSAA502 is a multi-purpose adapter that is used for eliminating ambient light or to provide a stable enhanced field of view at optimal working distance. Great for various applications such as checking the skin for dermatology or looking at a dogs ear for veterinarians. Simply place the Dino-Lite on top of the adapter and your ready to use.

Approximate field of view (FOV) and magnification using adapter:

Long working distance models: 1.5cm x 1.0cm at 35X

200x capable models: 4cm x 3cm at 10X

AM413TL-M40: 3cm x 2cm at 15X

