



## Model : SL-ZW1-Side Light with Polarizer

- Best results with AM413ZT, AM413ZTL or the AD413ZT.
- A polarized USB powered LED light
- Simple polarizer adjustment
- Built in illumination powered via USB

| Product Description |        |             |                           |
|---------------------|--------|-------------|---------------------------|
| Product Number      | SL-ZW1 | Dimension   | (L) 40cm x (W) 1.6cm      |
| Warranty Period     | 1      | Type        | Side Light with Polarizer |
| Condition           | New    | Unit Weight | 56g                       |

## Product Applications

The SL-ZW1 is not your ordinary articulate USB powered light. The SL-ZW1 is a polarized USB powered LED light that can be used with polarized Dino-Lite models.

With a simple rotation, you can adjust the amount of polarized illumination required for the right image.

To take advantage of the polarizing light, the accessory is best to be used with Dino-Lite models that have the polarization ability such as the AM413ZT, AM413ZTL or the AD413ZT.

This accessory is great for applications such as jewelry inspection, and also great for conditions where the Dino-Lite is too far from the object but still requires polarized illumination. If the Dino-lite does not have the polarization feature, you can still use it as additional external lighting or as a side light to enhance the view of surface topography! This desirable tool is a nifty accessory that is a great addition to the Dino-Lite digital handheld microscope family.

